



# Prototype Final P&ID



## Main Learnings

- Control  $\rightarrow$  high pressure & on/off solenoids is hard to make smooth & consistent when balancing IP, PEEP & B.P.M
- Flow & FiO<sub>2</sub> can be controlled through 2x flow meters on the air & oxygen lines
- System needs to be sealed well so it can control IP & PEEP better
- ~~System would operate much smoother  $\rightarrow$  variable flow control over a wide range; rather than pulsating on/off valves.~~
- Flow rate  $\rightarrow$  1-3 bar & irrigation solenoids is too high. need better flow control to reduce flow and make smooth control more manageable.

## Control Functions

- IP
- PEEP
- B.P.M
- ~~• FiO<sub>2</sub> (not really)~~
- 


## Measurement Functions

- TV
- Flow
- Pressure curve
- B.R.M & total breaths
- Minute ventilation
- RH (if humidifier used)
- Air temp

B.P.M: controlled using ~~either~~ flow rate & or  
time held at IP & PEEP.


IP: controlled using on/off value or variable flow rate

PEEP: controlled using on/off value or variable flow rate



- Graph above shows a possible way to smoothly control IP, PEEP & B.P.M with on/off solenoid values and some manual flow control. (See P&ID on next page).

## P&ID w Manual Flow Control



### Control logic

IP: SV-001 & 002 are opened if pressure is below IP. SV-001 & 002 close when IP is reached. SV-003 is closed throughout.

PEEP: SV-001 & 002 are closed. SV-003 is open if pressure is above PEEP. SV-003 is closed when pressure reaches PEEP.

B.P.M: This is equal to  $(\text{time taken to reach IP}) + (\text{time held at IP}) + (\text{time taken to reach PEEP}) + (\text{time held at PEEP})$ . The time held at IP & time held at PEEP can be controlled by the timing of SV-001, 002 & 003. The time to reach IP & PEEP can be manually set using OF-001, 002 & 003. Therefore B.P.M can be controlled unless  $(\text{time held at IP}) + (\text{time held at PEEP})$  needs to be negative.

I:E ratio: This ratio is equal to  $(\text{time taken to reach IP}) + (\text{time held at IP}) : (\text{time taken to reach PEEP}) + (\text{time held at PEEP})$ . This ratio can be adjusted by changing  $(\text{time held at IP}) : (\text{time held at PEEP})$  while maintaining B.P.M.

$F_iO_2$ : The oxygen concentration can be adjusted by altering the opening times of SV-001 & SV-002 ~~&~~ to ensure that ~~that~~ the flow rate ratio ensures the correct  $O_2$  concentration. See table below:  
 (assuming air has 20%  $O_2$  & Oxygen is 100%  $O_2$ )

~~(SV-001)~~ ~~(SV-002)~~  $x = \text{required } O_2 \text{ concentration}$

$$\textcircled{1} \quad QT-001 = \left( \frac{1-x}{0.8} \right)$$

$$\textcircled{2} \quad QT-002 = \left( \frac{1}{0.8} \left( 1 - \frac{1-x}{0.8} \right) \right)$$

Therefore for the following  $O_2$  concentrations the corresponding flow rate ratios will need to be obtained:

$O_2\% (x)$	$\frac{QT-002}{QT-001 + QT-002}$
20%	0
30%	0.125
40%	0.25
50%	0.375
60%	0.500
70%	0.625
80%	0.750
90%	0.875
100%	1

The  $O_2$  concentration would need to be maintained while ensuring that B.P.M, IP & PEEP are maintained.

## Interlocks

1. When SV-001 or SV-002 are open, SV-003 must be closed.
2. When SV-003 is open, SV-001 & SV-002 must be closed.
3. When SV-001 is open, SV-002 must be closed.
4. When SV-002 is open, SV-001 must be closed.
5. ~~QT-001 + QT-002 = QT-003~~ (for certain IP & PEEP take precedence over B.P.M. i.e. B.P.M can be automatically increased to ensure IP & PEEP are met. (Alarm should be raised)

## Alarms

SP alarms: If any items such as IP ~~less~~, PEEP, B.P.M & FiO<sub>2</sub> ever deviate from their set point values by more than specified amounts over a specified time frame then an alarm for that specific item should be raised.

Over pressure: If pressure in either PT-001 or PT-002 is recorded above a specified value (i.e. 60 cm H<sub>2</sub>O) then an alarm should be raised and ventilator shut down.

High B.P.M: If B.P.M goes above a set <sup>high</sup> value then an alarm should be raised.

Low B.P.M: If B.P.M goes below a set <sup>low</sup> value then an alarm should be raised.

TV asymmetric: Volume inhaled is ~~asymetric~~ different to volume exhaled (outside a certain margin). This can only be implemented if QT-003 is used.

TV + or - : Tidal volume increases or decreases (~~be~~ outside a certain margin) without any changes in other parameters.

## Patient triggered breathing

- If PT-001 (or PT-002) see a pressure <sup>a set amount</sup> lower than that of the PEEP set point. ~~with a margin~~ This should trigger an inspiration. ~~also~~ During patient triggered breathing <sup>a</sup> ~~the~~ set minimum B.P.M should still be met and an alarm raised if ~~the~~ a maximum B.P.M is exceeded. A minimum B.P.M will be ensured by having an underlying cycle running which will breath for the patient if they take too long to trigger.